

# Pulsed Xenon light source DB spectrophotometer



Windows Graphic Interface



News Release NL348

January 2020

Our new catalogue for our latest range of UV Vis instruments is now available plus dealer net price list

## The all new S-10 Superspec series UV - VIS - NIR Spectrophotometers

A superb new range of spectrophotometers providing low stray light, high sensitivity, extremely reliable and the best value on the market today

This is a great range of instruments economically priced, providing very high reliability and superb all round performance

More information at our website

[sales@spectrolab.co.uk](mailto:sales@spectrolab.co.uk)

## Superspec S-10

With Xenon light source One lamp all wavelengths

### Features

Double beam Spectrophotometer

CFR 21 compliant

PMT detector

Wavelength 190 – 1100nm

SBW Resolution 0.1nm

Working modes T.A.C.E.

Delivery 7 Days

[sales@spectrolab.co.uk](mailto:sales@spectrolab.co.uk)

