

News Letter 441 June 2022

The finest Atomic Absorption systems in the world

The Double beam AAS S44 Series from Spectrolab Atomic Absorption at its best

Spectrolab sell and support GBC products in many countries. The SavantAA is the 5th generation AAS It draws on a long tradition of excellence in design. The SavantAA atomic absorption spectrometer is a fully automated analyser capable of handling hundreds of samples per day. There



are three SavantAA models, the SavantAA, SavantAA Σ AAS and the SavantAA Z enduro. Powerful software allows unprecedented instrument control and data manipulation whilst following all the regulatory QC protocols and CFR 21 compliance. All instruments have USB communications for today's computer requirements.

For all models there is a complete range of accessories available to further expand the capability of these instruments including graphite furnace system, hydride system, mercury concentrator, HCL lamps and auto samplers, auto diluters etc.

All controlled by fully featured Windows[®] based software.

And the price ? From just \$8000

Best value on the market and best technology



sales@spectrolab.eu www.spectrolab.co.uk