



SPECTROLAB

Berkshire England RG17 0NL UK C 00362

sales@spectrolab.eu sales@spectrolab.co.uk www.spectrolab.co.uk

News letter 453 August 2022

MWD-6100 Series Microwave digestion systems

MWD-6100 series is our latest series of smart microwave digesters designed for both high performance and safe operator protection. The digestors are perfect for ICP and many other applications. They are easy to use and offer a very high level of operator safety



Features

- Optional vessel quantity, 6 to 12 vessels to meet different digestion requirements.
- Peek and Glass construction
- Real-time monitoring of both temperature and pressure in each vessel.
- Contactless sensor monitoring with no burst disk design saves consumables cost and maintenance cost.
- 7-inch colour Touch screen with user-friendly interface for easy operation.
- CFRP outer shell material with high strength ensures impact resistance.
- Pre-installed general standard methods, users can also create, save, modify and delete the method.
- 316 L industrial stainless-steel tank with multilayer Teflon coating avoids acid corrosion, also

improved cooling efficiency.

Best value on the market today

Contact us at

sales@spectrolab.co.uk